


Olis GlobalWorks 3D Analysis Software

Post-collection software with a modern emphasis on multiple wavelength -- Global -- data analysis and handling.



Global analysis works better, providing more accurate rate constants and definitive spectral reconstruction.

Single wavelength analysis is often incomplete and incorrect. But analyzing large **multiple wavelength** data sets is beyond the scope of most software packages or too slow to be practical.

In the early 1990s, after introducing a rapid-scanning spectrophotometer which produces 1 MB of data per second, Olis programmers were charged with developing fast and effective **multiple wavelength** data analysis software. Today's Olis GlobalWorks software is the extraordinary result of more than 14 years of implementing and improving algorithms for 3D analysis.


Start with **multiple wavelength** data (spectra) acquired as a function of time, temperature, or other process.


Apply SVD. Instantly, the spectral species are found and separated from the noise in the raw data.


Now, using precise graphic and numeric parameters returned by SVD to conclude how many species are present, choose a chemical model to define the reaction (e.g., $A \Leftrightarrow B \to C$). Moments elapse from raw 3D data entry to fitted results.


Do not let insufficient data or frustrating software limit you to single wavelength data analysis. Move to Olis GlobalWorks. More than 70 kinetic and equilibrium cases are supported for files containing ten to thousands of scans.

Step 1: Time, Temperature, or Other Process Dependent 3D Data













Step 3: Chemical Models for Chosen Number of Species








Answer: Numeric Results, and plotted Residuals, Kinetics and Spectra.




2D is Insufficient

Seminal example¹ of 'simple' two rate reaction requiring multiple wavelength data for correct analysis.



Synthetic 3D Data

A multiple wavelength data set containing two species, A and B, with associated rate constants of 5 and 6, was synthesized; a modest amount of noise was added.




Fixed Wavelength Trace

Looking only at the kinetic trace at a given wavelength, there is no clear suggestion of two rates. This looks like a single exponential. And any one wavelength analyzed is insufficient and even misleading, as we see on the facing page.

¹Knutson JR., 1992 Methods Enzymol. 210 p 357-74. Maeder, M., and Zuberbuhler, A.D., 1990 Anal. Chem. 62 p 2220.

😽 GlobalWorks - [Data Display] 📆 File Edit View Fits Tools Chart Help E I BRAAC X (X 14 □ | 4 4 17 € | SM 0.000 (1.3066±0.0008)*E+00 (5.551±0.013)*E+00 -0.010 0.00.10.20.30.40.50.60.70.80.91.0 Time, sec 1.0 Absorbance 0.4 0.1 0.7 1.0 0.0 0.2 0.3 0.4 0.5 0.6 0.8 0.9 Time, sec



Fit at 600 nm

When one considers data extracted at 600 nm - a completely reasonable wavelength to choose - the single rate returned is 5 s⁻¹. The residuals give no clue that the fit is more complicated than we show.

Fit at 500 nm

When one considers data extracted at 500 nm - again, a completely reasonable wavelength to choose - the single rate returned is 5.5 s⁻¹. Again, with no other information, we would confidently declare this is one exponential.


Fit at 450 nm

This time, your alert should go up (but would it?) because the rate constant is not within the 5–5.5 s¹ range that the other fits calculated, but is now at 6 s¹! What is wrong? The fit in all three cases - if considered individually - would seem conclusive. The rate is 5 s¹. Or 5.5 s¹. Or 6 s¹. But when they are put together, and if a knowledgeable kineticist is interpreting the results, the conclusion would have to be made that "More experiments need to be done!" But, oops. The enzyme is all gone. Or, the student has graduated. Or...?

GlobalWorks - [Data Display]


3D is Required & Successful

Continuing with the previous example, see how 3D global analysis easily finds two rates with rate contants of 5 and 6.



Synthetic 3D file

The synthetic data set was created to have 1000 multiple wavelength data points, exactly as would be produced by an Olis RSM 1000 rapid-scanning spectrophotometer (which collects 1,000 scans per second).




Eigenvector display

Hand all or part of the raw data to the OLIS Global Fitting scheme. The first step is singular value decomposition (SVD), or 'factor analysis.' Factor analysis finds things with some shape ("spectra") which are changing at some rate. The results are shown in "eigenvectors," graphical representations of shapes. The upper row shows spectral eigenvectors; the lower row shows kinetic eigenvectors. Eigenvectors without structure are "noise." Thus, noise contributions are isolated from the useful information and can be excluded from subsequent data handling.


Possible Chemical Equations

With the proposed number of species selected, one moves to defining how these species are related. Do we have one form of A changing into another form (as happens in a protein-fold reaction)? Or do we have a bleaching reaction? Are the species related sequentially (A → B)? These are questions we cannot answer without spectra.




The Answer: 2 Species, 2 Rates

The data were fitted to selection D: (A →..., B →) Since these are synthetic data, we know the answer. The known rates of 5.0 and 6.0 are found exactly by Olis GlobalWorks.




3D Data Supports Determination of the Order of Rate Constants

With single wavelength data, the order of two rate constants in a simple sequential mechanism cannot be determined; with multiple wavelength data, it can!



Raw 3D Data file

The synthetic data set was created to have 1000 multiple wavelength data points. Rate constants of 5 sec¹ and 7 sec¹ were used for this $A \rightarrow B \rightarrow C$ case.



Results of SVD

Eigenvectors suggest (clearly) three species undergoing changes, and noise.

Chemical models


Fit options available when there are three species and two or three rates.

	Select Fit Model	
	Sequential Cases	
A-> B-> C Fast/Slow Rate	3 Species Sequential	
A-> B-> C Slow/Fast Rate	3 Species Sequential	
A + A -> B -> C	3 Species Sequential, First Process Second Order	
A-> B, B + B -> C	3 Species Sequential, Second Process Second Order	
	Multiple First Order Cases	
A ->, B -> C	First Order Decay and Independent First Order Growing in	
A -> B, C ->	First Order in Growing and Independent First Order Decay	
A -> B, C -> B	2 Independent First Orders, Same Product	
A ->, B->, C	2 Independent First Orders and Background	
A →, B →, C →	->, B ->, C -> 3 Independent First orders	
> A,> B, C 2 First Orders Growing in and Background		
	Second Order Cases	
A + B -> C	Heterogeneous Second Order	
A -> B, B + B -> C + A	First Order Growing in followed by Second Order	
A -> B, C ->, B + C -> A +	Mixed Reduction Reactions	
	Fit Models Without Unique Solutions	

Fit Model A-> B-> C, Fast/Slow Rate Constants

Incorrect Fit


Since we synthesized the data, we know that there should not be a shoulder at 400 nm in spectrum B. Thus, the fast rate cannot be first.




Fit Model A-> B-> C Slow/Fast Rate Constants


Correct Fit

While the kinetic results are identical with this fit and the previous (notice the residuals and the kinetic traces), the spectra are now correct with the slow rate associated with k1, not k2.







Apply SVD



Choose Fit From One Species Mechanisms



Evaluate Results



Proving the Absolute Absence of Bias in GlobalWorks

Using the Olis RSM 1000 in its fluorescence mode, 50 scans were collected in 0.05 seconds from light emitted by overhead fluorescent lamps. The mercury lines and the 120 Hz modulation of the lamp jump from in the eigenvectors (second panel) and the answer (final graphs).

SVD identified the changes in the light caused by the lamp's on/off operation at 60 Hz, and the spectrum of the light (a continuum with sharp Hg emission lines). SVD acts with no bias about the source of the data, chemical, kinetic, thermal, or other process. And, since user starting values are only used for the most complicated cases, bias cannot be entered during the fit, either.

Collection: 50 scans collected in 0.05 seconds of light emitted by overhead fluorescent lamps

SVD: Single species

Kinetic Fit: Select 1 species and then fit single rate so as to construct the spectrum

Mechanism: Not a kinetic process

Comments: Fluorescent lights switching on and off at 120 Hz, demonstration of a particular shape (the spectrum) varying in a particular way (sinusoidal). The eigenvectors are presented in the way shown to emphasize that the kinetic and spectral eigenvectors are related. One reads the display by noting that the kinetic eigenvectors show a particular time course and that it is the corresponding spectral eigenvector which varies in that way.

These data, as with all data shown in this brochure, are provided for pedagogical, demonstration, and testing purposes in the Olis GlobalWorks software.

Temperature Dependent Data

Using an Olis DSM CD spectrophotometer with automatic temperature regulation, 38 CD scans were collected over 50 degrees.

SVD suggests two species, which, are one form of a protein turning into another form of the protein in a (potentially) reversible manner.

As a thermal study, numeric results returned include the transition ("melting") temperature and the enthalpy value(s).

With single wavelength results only, one could say nothing about the substance undergoing the thermal denaturation. With multiple wavelength results, one has the correct spectra of the native and unfolded form of the human serum albumin in addition to the more accurate determination of the transition temperature. How much more satisfying this full result is!


Collection: 38 CD scans collected as a function of temperature with an OLIS DSM CD spectrometer.


SVD: Shows clearly that two species are involved. Note that the X- axis is temperature, not time.

Thermodynamic Fit: Calculate thermodynamic properties (e.g. melting temperature and enthalpy).


Mechanism: N↔U

Comments: These data show that the protein changed as the temperature was increased. A fit to these data will give a value for the 'melting temperature' or transition temperature. Spectral reconstruction produced the spectrum of the native (folded) and denatured (unfolded or melted) forms.






Choose Fit From Equilibrium Mechanisms



Evaluate Results



These data, as with all data shown in this brochure, are provided for pedagogical, demonstration, and testing purposes in the Olis GlobalWorks software.

Currently Supported Mechanisms


Equilibrium Mechanisms				
	Mative to Unfolded State	Constant Enthalpy with Temperature		
N→U		Constant Enthalpy with Temperature		
N→U	Native to Unfolded State	Enthalpy Varies with Temperature CD Titration Curve as function of Added Volume		
N→U CD=Cdo[V(V+Vo)] Returns Vo				
N↔U	Native to Unfolded State	Constant Enthalpy with Temperature		
N↔U Native to Unfolded State		Enthalpy Varies with Temperature		
N.N↔U Dimer Native to Unfolded State		Constant Enthalpy with Temperature		
N.N↔U	Dimer Native to Unfolded State	Enthalpy Varies with Temperature		
N.N.N↔U	Trimer Native to Unfolded State	Constant Enthalpy with Temperature		
N⇔I↔U	Native to Unfolded State with Intermediate	Constant Enthalpy with Temperature		
N⇔l1⇔l2⇔U	Native to Unfolded State with 2 Intermediates	Constant Enthalpy with Temperature		
CD=CDo*[V/(V+Vo)]	Returns Vo	CD Titration Curve as Function of Added Volume		
F = 1/(1+10^(pH-pK))		Henderson-Hasselbach Eqn. to Determine pK		
N+D↔U	Denatruration Process	Eqbm. Const. Ko Determiined, N Constant, D Varied		
Spectral Changes as a		Delta V, Eqbm. Const. Ko Determined, T Constant		
One Species Mechanisms				
	Sequenti			
A→		Simple First Order Decay		
A + A→		Second Order Decay		
→A		Growing First In		
Two Species Mechani				
	Sequenti			
A→B		First Order with Growing In		
A↔B		Reversible Reaction		
Multiple First Order Cases				
A→, B		First Order with Background		
A→, B→		2 First Order Decays		
A→, B→	Linear Drift	First and Zeroth Orders		
→A, B		First Order Growing In with Background		
→A, B→		First Order Growing In and Decay		
Second Order Cases				
A+A→B		Second Order with Growing In		
A+A→B		Second Order Decay with Background		
A→, A+A→B		Competing First and Second Orders with Background		
A→B, A+A→B		Competing First and Second Orders with Same Product		
Three Species Mechanisms				
Sequential Cases				
A→B→C	Fast/Slow Rate Constants	3 Species Sequential		
A→B→C	Slow/Fast Rate Constants	Special 3 Species Sequential		
A+A→B→C		3 Species Sequential, First Process Second Order		
A→B, B+B→C		3 Species Sequential, Second Process Second Order		
	Multiple First	Order Cases		
A→, B→C		First Order Decay & Independent First Order Growing In		
A→B, C→		First Order & Independent First Order Decay		
A→B, C→B		2 Independent First Orders, Same Product		
A→, B→, C		2 Independent First Orders and Background		
A→, B→, C→		3 Independent Fisrt Orders		
→A,→B, C 2 First Orders Growing In and Background				

Currently Supported Mechanisms

Second Order Cases				
A+B→C	Heterogeneous Second Order			
A→B, B+B→C+A	First Order In followed by Second Order			
A→B, C→, B+C→A	Mixed Reduction Reactions			
Fit Models Without Unique Solutions				
A → B → C Partially Reversible Rise-Fall				
A→B↔C	Partially Reversible Rise-Fall			
Four Species Mechanisms				
Sequential Cases				
$A\rightarrow B\rightarrow C\rightarrow C$	4 Species Sequential			
Multiple First Order Cases				
A→, B→C	First Order Decay & Independent First Order Growing In			
A→B, C→	First Order In Growing & Independent First Order Decay			
Second Order and Enzyme Cases				
A→, B→, C→, D	Pulse Radiolysis Case			
E+S→ES→E+P	Irreversible Enzyme Case			
Fit Models Without Unique Solutions: Reversible Sequential Cases				
$A \leftrightarrow B \rightarrow C \rightarrow D$	Sequential with Step 1 Reversible			
A⇔B⇔C⇔D	Fully Reversible Sequential			
A↔B, B↔C, B↔D	Branching Reaction			
Fit Models Without Unique S	Solutions: Reversible Enzyme Cases			
E+S↔ES→E+P	Partially Reversible Enzyme Reaction			
E+S+ES+P	Fully Reversible Enzyme Reaction			
Five Species Mechanisms				
Sequential Cases				
$A\rightarrow B\rightarrow C\rightarrow D\rightarrow E$	5 Species Sequential			
Multiple I	First Order Cases			
A→, B→, C→, D→, E	4 Independent First Orders and Background			
A→, B→, C→, D→, E→	First Order In Growing & Independent First Order Decay			
Enzyme Case				
E+S→ES→EP→E+P	Enzyme Reaction, 2 Enzyme Complexes			
Fit Models Without Unique Solutions: Reversible Sequential Cases				
A↔B↔C↔D↔E	Fully Reversible Sequential			
	Fully Reversible Sequential Solutions: Reversible Enzyme Cases			
Fit Models Without Unique S	Solutions: Reversible Enzyme Cases			
Fit Models Without Unique S E+S↔ES→EP→E+P	Partially Reversible Enzyme Reaction, 2 Enzyme Complexes			
Fit Models Without Unique S E+S↔ES→EP→E+P E+S↔ES↔EP↔E+P Six Species Mechanisms	Partially Reversible Enzyme Reaction, 2 Enzyme Complexes			
Fit Models Without Unique S E+S↔ES→EP→E+P E+S↔ES↔EP↔E+P Six Species Mechanisms	Partially Reversible Enzyme Reaction, 2 Enzyme Complexes Fully Reversible Enzyme Reaction, 2 Enzyme Complexes			
Fit Models Without Unique S E+S↔ES→EP→E+P E+S↔ES↔EP↔E+P Six Species Mechanisms Sequ A→B→C→D→E→F	Partially Reversible Enzyme Reaction, 2 Enzyme Complexes Fully Reversible Enzyme Reaction, 2 Enzyme Complexes Fully Reversible Enzyme Reaction, 2 Enzyme Complexes			
Fit Models Without Unique S E+S↔ES→EP→E+P E+S↔ES↔EP↔E+P Six Species Mechanisms Sequ A→B→C→D→E→F Multiple I A→, B→, C→, D→, E→, F	Partially Reversible Enzyme Reaction, 2 Enzyme Complexes Fully Reversible Enzyme Reaction, 2 Enzyme Complexes Fully Reversible Enzyme Reaction, 2 Enzyme Complexes Itential Cases 6 Species Sequential			
Fit Models Without Unique S E+S↔ES→EP→E+P E+S↔ES↔EP↔E+P Six Species Mechanisms Sequ A→B→C→D→E→F Multiple I	Partially Reversible Enzyme Reaction, 2 Enzyme Complexes Fully Reversible Enzyme Reaction, 2 Enzyme Complexes Fully Reversible Enzyme Reaction, 2 Enzyme Complexes Tential Cases 6 Species Sequential First Order Cases			
Fit Models Without Unique S $E+S \leftrightarrow ES \rightarrow EP \rightarrow E+P$ $E+S \leftrightarrow ES \leftrightarrow EP \leftrightarrow E+P$ $Six Species Mechanisms$ $Sequential A \rightarrow B \rightarrow C \rightarrow D \rightarrow E \rightarrow F$ $Multiple B \rightarrow B $	Partially Reversible Enzyme Reaction, 2 Enzyme Complexes Fully Reversible Enzyme Reaction, 2 Enzyme Complexes Fully Reversible Enzyme Reaction, 2 Enzyme Complexes Idential Cases 6 Species Sequential First Order Cases 5 Independent First Orders and Background			

This supported mechanism list continues to grow as people ask for additional equations. If you do not see the mechanism which may apply to your multiple wavelength data, provide us the equation and we will incorporate it in Olis GlobalWorks at no cost to you. If you can express your equation algebraically, we can create a fit for it.

Numeric Results Returned



Quintessentially New Global Fitting Software Uses Simplex Method and Matrix Exponentiation

Analyzing data in a multidimensional form is superior to analyzing 2D data. The new GlobalWorks algorithms by Olis make analyzing multidimensional data easy and fast. Algorithms for data as a function of time, temperature, pressure, concentration, and angle are now supported; additional models can be added.

Matrix Exponential can be extended to very complicated cases, including reversible cases for which solutions are algebraically intractable. The limit to fit complexity is not computational, but the data (the more rates to be determined, the higher the S/N of the data must be).

There is no longer justification to limit one's analyses to single wavelength data. Use multiple wavelength data and Olis GlobalWorks for better rate constants, digital noise reduction, and chemical mechanism testing. Results are so much more satisfying and defensible!

Multidimensional data provide a much stronger basis for accurate conclusions.

GlobalWorks provides
the best algorithms
extant for
instantaneous and
correct analysis of
multidimensional data.

Correct Results the Multidimensional Way

Olis GlobalWorks differs from all other kinetic fitting software in three ways. Our SVD is 1300 fold faster than the original algorithm. Instead of Levenberg-Marquardt, we use Downhill Simplex. And instead of numeric integration routines such as Runge-Kutta, we use Matrix Exponentiation.

Matrix Exponential can be extended to very complicated cases, including reversible cases for which solutions are algebraically intractable. Consider the simple case of $A \rightarrow B$, such as defines a protein unfold experiment. By applying the Law of Mass Action, a matrix containing only rate constants and zeros is produced. The columns represent the concentrations and the rows their derivatives with respect to time. The matrix is exponentiated to provide the solution for $A \rightarrow B$.

For additional details on how Olis GlobalWorks makes 3D analysis fast, easy, and indispensable, see Methods in Enzymology, volume 384, 2004, chapters 1, 2, and 3. The three data sets used in this document (protein unfolding, Xylenol reacting with iron, and nucleic acids undergoing thermal denaturation) are some of the demonstration files provided with the software.

Use the software today. A 30 day free trial version is available at http://www.olisweb.com/software/

For more information on this and other Olis products:

Visit www.olisweb.com

Write sales@olisweb.com

Call 1-800-852-3504 in the US & Canada

1-706-353-6547 worldwide