OLIS SLM 8000

Highly esteemed benchmark spectrofluorometer

Our modernization procedure will add decades of reliable operations from your long used and trusted SLM. All hardware is replaced with modern. The ageless detectors are retained with all new circuitry. Original accessories can be retained, computerized, or replaced as needed. New accessories can be added, giving you the exact instrument you need. And, of course, all software is the OLIS Windows 10 SpectralWorks, a comprehensive data acquisition and instrument control program.



<u>Standard Acquisition Mode:</u> Fluorescence

Enhancements Supported:

Peltier Thermal Control Stopped Flow Titrator

Process for Upcycling

- 1. We ship a crate to you that you'll pack the SLM in
- 2. After we receive the instrument, we replace 100% of the electronics with the OLIS Master Control Box and Win 10 computerization
- 3. We exhaustively test the unit for performance with the accessories you have or chose to purchase from us; and ship the unit back
- 4. Alternatively, you can receive the upgrade kit from us and do the simple replacement with our virtual support.