



OLIS Software

Three programs: GlobalWorks, both freestanding and bundled with data collection module; OlisWorks, data collection module for the HP 8453 and CLARITY CCD; and SpectralWorks, the data collection module for most OLIS spectrophotometers.




OlisWorks

A major advance over ChemStation for the HP 8453 diode array. Also the controlling interface for our only CCD spectrometer, the CLARITY CCD.



SpectralWorks

Comprehensive data acquisition and instrument control program for 2D and 3D data files collection as a function of time, temperature, titration, pressure, concentration, or other variables on all scanning OLIS spectrophotometers.



GlobalWorks

Post-collection module with a modern emphasis on multiple wavelength – global – data analysis and handling. SVD for 3D data fitting and spectral reconstruction, plus dozens of algorithms for 2D data.

